


Product Name: Recombinant Human MAP2K4 Protein
Catalog #: PEH7261



Summary

Name	MP2K4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P45985
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAPSPSGGG GSGGGSGSGT PGPVGSPAPG HPAVSSMQGK RKALKLNFAN PPFKSTARFT LNPNTGVQN PHIERLRTHS IESSGKLKIS PEQHWDFTAE DLKDLGEIGR GAYGSVNKMV HKPSGQIMAV KRIRSTVDEK EQKQLMDLD VVMRSSDCPY IVQFYGALFR EGDCWICMEL MSTSFDFYK YVYSVLDDVI PEEILGKITL ATVKALNHLK ENLKIIHRDI KPSNILLDRS GNIKLCDFGI SGQLVDSIAK TRDAGCRPYM APERIDPSAS RQGYDVRSDD WSLGITLYEL ATGRFPYPKW NSVFDQLTQV VKGDPPQLSN SEEREFSPSF INFVNLCLTK DESKRPKYKE LLKHPFILMY EERAVEVACY VCKILDQMPA TPSSPMYVD

image



Background

Alternative Names	Dual specificity mitogen-activated protein kinase kinase 4 (MAP kinase
--------------------------	--

Product Name: Recombinant Human MAP2K4 Protein
Catalog #: PEH7261



kinase 4) (MAPKK 4) (EC 2.7.12.2) (JNK-activating kinase 1) (MAPK/ERK kinase 4) (MEK 4) (SAPK/ERK kinase 1) (SEK1) (Stress-activated protein kinase kinase 1) (SAPK kinase 1) (SAPKK-1) (SAPKK1) (c-Jun N-terminal kinase kinase 1) (JNKK)

Note

For Research Use Only , Not for Diagnostic Use