


Product Name: Recombinant Human KCNJ3 Protein
Catalog #: PEH7265



Summary

Name	KCNJ3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48549
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERHNSVECLD GLDDITTKLP SKLQKITGRE DFPKKLLRMS STTSEKAYSL GDLPMKLQRI SSVPGNSEEK LVSKTTKMLS DPMSQSVADL

image



Background

Alternative Names	G protein-activated inward rectifier potassium channel 1 (GIRK-1) (Inward rectifier K(+) channel Kir3.1) (Potassium channel, inwardly rectifying subfamily J member 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use