


**Product Name: Recombinant Human CSNK2A1 Protein**  
**Catalog #: PEH7282**



## Summary

<b>Name</b>	CSK21_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P68400
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VEWGNQDDYQ LVRKLGRGKY SEVFAINIT NNEKVVVKIL KPVKKKKIKR EIKILENLRG GPNII TLADI VKDPVSRTPA LVFEHVNNTD FKQLYQTLTD YDIRFYMYEI LKALDYCHSM GIMHRDVKPH NVMIDHEHRK LRLIDWGLAE FYHPGQEYNV RVASRYFKGP ELLVDYQMYD YSLDMWSLGC MLASMIFRKE PFFHGHNDYD QLVRIAKVLG TEDLYDYIDK YNIELDPRFN DILGRHSRKR WERFVHSENQ HLVSPEALDF LDKLLRYDHQ SRLTAREAME HPYFYTVVKD QARMGSSSMP

## image



## Background

**Alternative Names** Casein kinase II subunit alpha (CK II alpha) (EC 2.7.11.1)

**Product Name: Recombinant Human CSNK2A1 Protein**  
**Catalog #: PEH7282**



---

**Note**

For Research Use Only , Not for Diagnostic Use