


Product Name: Recombinant Human NFkB2 Protein
Catalog #: PEH7286



Summary

Name	NFKB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q00653
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VIVEQPKQRG FRFRYGCEGP SHGGLPGASS EKGRKTYPTV KICNYEGPAK IEVDLVTHSD PPRAHAHSLV GKQCSELGIC AVSVGPKDMT AQFNNLGV LH VTKKNMMGMT IQKLQRQLR SRPQGLTEAE QRELEQEAKE LKKVMDLSIV RLRFSAFLRA SDGSFSLPLK PVISQPIHDS KSPGASNLKI SRMDKTAGSV RGGDEVYLLC DKVQKDDIEV RFYEDDENGW QAFGDFSPTD VHKQYAIVFR TPPYHKMKIE RPVTVFLQLK RKRGGDVSDS KQFTYYPLVE DKEEVQRKRR

image



Background

Alternative Names	Nuclear factor NF-kappa-B p100 subunit (DNA-binding factor KBF2) (H2TF1) (Lymphocyte translocation chromosome 10 protein) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2) (Oncogene Lyt-10) (Lyt10)
--------------------------	--

Product Name: Recombinant Human NFkB2 Protein
Catalog #: PEH7286



[Cleaved into: Nuclear factor NF-kappa-B p52 subunit]

Note

For Research Use Only , Not for Diagnostic Use