


## Summary

<b>Name</b>	ANK2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q01484
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	35 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DSGEKFNQSS QRRKRPKKSD SNASFLRAAR AGNLDKVVVEY LKGGIDINTC NQNGLNALHL AAKEGHVGLV QELLGRGSSV DSATKKGNTA LHIASLAGQA EVLKVLVKEG ANINAQSQNG FTPLYMAAQE NHIDVVKYLK ENGANQSTAT EDGFTPLAVA LQQGHNQAVA ILLENDTKGK VRLPALHIAA RKDDTKSAAL LLQNDHNADV QSKMMVNRTT ESGFTPLHIA AHYGNVNVAT LLNRGAADV FTARNGITPL HVASKRGNTN MVKLLDRGG QIDAKTRDGL TPLHCAARSG

## image



## Background

**Alternative Names** Ankyrin-2 (ANK-2) (Ankyrin-B) (Brain ankyrin) (Non-erythroid ankyrin)

**Product Name: Recombinant Human ANK2 Protein**  
**Catalog #: PEH7287**

---



**Note**

For Research Use Only , Not for Diagnostic Use