

Product Name: Recombinant Human CACNA1D Protein
Catalog #: PEH7288



Summary

Name	CAC1D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01668
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTEGNLEQAN EELRAVIKKI WKKTSMKLLD QVVPPAGDDE VTVGKFYATF LIQDYFRKFK KRKEQGLVGK YPAKNTTIAL QAGLRTLHDI GPEIRRAISC DLQDDEPEET KREEEDDVFK RNGALLGNHV NHVNSDRRDS LQQTNTTHRP LHVQRPSIPP ASDTEKPLFP PAGNSVCHNH HNHNSIGKQV PTSTNANLNN ANMSKAAHGK RPSIGNLEHV

image



Background

Alternative Names	Voltage-dependent L-type calcium channel subunit alpha-1D (Calcium channel, L type, alpha-1 polypeptide, isoform 2) (Voltage-gated calcium channel subunit alpha Cav1.3)
--------------------------	--

Product Name: Recombinant Human CACNA1D Protein
Catalog #: PEH7288



Note

For Research Use Only , Not for Diagnostic Use