


Product Name: Recombinant Human CACNB1 Protein
Catalog #: PEH7290



Summary

Name	CACB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02641
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYPSPQEIPM EVFDPSQGK YSKRKGRFKR SDGSTSSDTT SNSFVRQGSA ESYTSRPSDS DVSLEEDREA LRKEAERQAL AQLEKAKTKP VAFAVRTNVG YNPSPGDEVP VQGVAITFEP KDFLHIKEY NNDWWIGRLV KEGCEVGFIP SPVKLDSLRL LQEQLRQNR LGSSKSGDNS SSSLGDVVTG TRRPTPPASA KQKQKSTEHV

image



Background

Alternative Names Voltage-dependent L-type calcium channel subunit beta-1 (CAB1) (Calcium channel voltage-dependent subunit beta 1)

Note

Product Name: Recombinant Human CACNB1 Protein
Catalog #: PEH7290



For Research Use Only , Not for Diagnostic Use