


Product Name: Recombinant Human BTK Protein
Catalog #: PEH7292



Summary

Name	BTK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q06187
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKRSQQKKKT SPLNFKKRLF LLTVHKLSYY EYDFERGRRG SKKGSIDVEK ITCVETVVPE KNPPPERQIP RRGEESSEME QISIIEFPY PFQVVYDEGP LYVFSPTTEL RKRWIHQLKN VIRYNSDLVQ KYHPCFWIDG QYLCCSQTAK NAMGCQILEN RNGSLKPGSS HRKTKKPLPP TPEEDQILKK

image



Background

Alternative Names	Tyrosine-protein kinase BTK (EC 2.7.10.2) (Agammaglobulinemia tyrosine kinase) (ATK) (B-cell progenitor kinase) (BPK) (Bruton tyrosine kinase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BTK Protein
Catalog #: PEH7292

