


Product Name: Recombinant Human BCL2L1 Protein
Catalog #: PEH7294



Summary

Name	B2CL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q07817
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFLSYKLSQK GYSWSQFSDV EENRTEAPEG TESEMETPSA INGNPSWHLA DSPAVNGATG HSSSLDAREV IPMAAVKQAL REAGDEFELR YRRAFSDLTS QLHITPGTAY QSFEQVVNEL FRDGVNWGRI VAFFSFGGAL CVESVDKEMQ VLVSRIAAMW ATYLNHLEP WIQENGWDT FVELYGNNA

image



Background

Alternative Names Bcl-2-like protein 1 (Bcl2-L-1) (Apoptosis regulator Bcl-X)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BCL2L1 Protein
Catalog #: PEH7294

