

Summary

Name	PDE4D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08499
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDVLAKELED VNKWGLHVFR IAELSGNRPL TVIMHTIFQE RDLLKTFKIP VDLITYLMT LEDHYHADVA YHNNIHAADV VQSTHVLLST PALEAVFTDL EILAAIFASA IHDVDHPGVS NQFLINTNSE LALMYNDSSV LENHHLAVGF KLLQEENCDI FQNLTKKQRQ SLRKMVIDIV LATDMSKHMN LLADLKT MVE TKKVTSSGVL LLDNYS DRIQ VLQNMVHCAD LSNPTKPLQL YRQWTD RIME EFRQGD RER ERGMEISPMC DKHNASVEKS QVGFIDYIVH PLWETWADLV HPDAQDILDT LEDNREWYQS

image



Background

Alternative Names	cAMP-specific 3',5'-cyclic phosphodiesterase 4D (EC 3.1.4.53) (DPDE3) (PDE43)
--------------------------	---

Product Name: Recombinant Human PDE4D Protein
Catalog #: PEH7297



Note

For Research Use Only , Not for Diagnostic Use