


Product Name: Recombinant Human DPYD Protein
Catalog #: PEH7299



Summary

Name	DPYD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q12882
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEQDETGKWN EDEDQMVHLK ADVVISAFGS VLSDPKVKEA LSPKFNRWG LPEVDPETMQ TSEAWVFAGG DVVGLANTTV ESVNDGKQAS WYIHKYVQSQ YGASVSAKPE LPLFYTPIDL VDISVEMAGL KFINPFGLAS ATPATSTSMI RRAFEAGWGF ALTKTFSLDK DIVTNVSPRI IRGTTSGPMY GPGQSSFLNI ELISEKTAAY WCQSVTELKA DFPDNIVIAS IMCSYNKNDW TELAKKSEDS GADALELNLS CPHGMGERGM GLACGQDPEL VRNICRWVRQ AVQIPFFAKL TPNVTDIVSI ARAAKEGGAN GVTATNTVSG

image



Background

Alternative Names	Dihydropyrimidine dehydrogenase [NADP(+)] (DHPDHase) (DPD) (EC 1.3.1.2) (Dihydrothymine dehydrogenase) (Dihydrouracil dehydrogenase)
--------------------------	---

Product Name: Recombinant Human DPYD Protein
Catalog #: PEH7299



Note

For Research Use Only , Not for Diagnostic Use