


Product Name: Recombinant Human PDE3B Protein
Catalog #: PEH7302



Summary

Name	PDE3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13370
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KESFKLMETQ QEEETKEDS RKLFEQGDW LEEAQSEQQ TNIEQEVSLD LILVEEYDSL IEKMSNWNFP IFELVEKMGE KSGRILSQVM YTLFQDTGLL EIFKIPTQQF MNYFRALENG YRDIPYHNRI HATDVLHAVW YLTTRPVPGL QQIHNGCGTG NETDSDGRIN HGRIAYISSK SCSNPDESYG CLSSNIPALE LMALYVAAAM HDYDHPGRTN AFLVATNAPQ AVLYNDRSVL ENHHAASAWN LYLSRPEYNF LLHLDHVEFK RFRFLVIEAI LATDLKKHFD FLAEFNAKAN DVNSNGIEWS NENDRLLVCQ VCIKLADING PAKVRDLHLK WTEGIVNEFY EQGDDEANLG LPISPFMDRS SPQLAKLQES FITHIVGPLC NSYDAAGLLP GQWLEAEEDN DTESGDDEDG EELDTEDDEM ENNLNPKPPR RKSRRRIFCQ LMHHLTENHK IWKEIVEEEE KCKADGNKLQ VENSSLPQAD EIQVIEEAD

image



Background

Product Name: Recombinant Human PDE3B Protein
Catalog #: PEH7302



Alternative Names

cGMP-inhibited 3',5'-cyclic phosphodiesterase B (EC 3.1.4.17) (CGIPDE1)
(CGIP1) (Cyclic GMP-inhibited phosphodiesterase B) (CGI-PDE B)

Note

For Research Use Only , Not for Diagnostic Use