


Product Name: Recombinant Human ROCK1 Protein
Catalog #: PEH7303



Summary

Name	ROCK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13464
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQYFELTQES KKAASRNQEQ ITDKDHTVSR LEEANSMLTK DIEILRRENE ELTEKMKKAE EYKLEKEEE ISNLKAAFEK NINTERTLKT QAVNKLAEIM NRKDFKIDRK KANTQDLRKK EKENRKLQLE LNQEREKFNQ MIVVKHQKELN

image



Background

Alternative Names	Rho-associated protein kinase 1 (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-35) (Rho-associated, coiled-coil-containing protein kinase 1) (Rho-associated, coiled-coil-containing protein kinase I) (ROCK-I) (p160 ROCK-1) (p160ROCK)
--------------------------	---

Note

Product Name: Recombinant Human ROCK1 Protein
Catalog #: PEH7303



For Research Use Only , Not for Diagnostic Use