


Product Name: Recombinant Human NFATC2 Protein
Catalog #: PEH7304



Summary

Name	NFAC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13469
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSSQSGSYEL RIEVQPKPHH RAHYETEGSR GAVKAPTGGH PVVQLHGYME NKPLGLQIFI GTADERILKP HAFYQVHRIT GKTVTTTSYE KIVGNTKVLE IPLEPKNNMR ATIDCAGILK LRNADIELRK GETDIGRKNT RVRLVFRVHI PESSGRIVSL QTASNPIECS QRSARELPMV ERQDTSCLV YGGQQMILTG QNFTSESKVV FTEKTTDGQQ IWEMEATVDK DKSQPMLFV EIPEYRNKHI RTPVKVNFYV INGKRKRSQP

image



Background

Alternative Names	Nuclear factor of activated T-cells, cytoplasmic 2 (NF-ATc2) (NFATc2) (NFAT pre-existing subunit) (NF-ATp) (T-cell transcription factor NFAT1)
--------------------------	--

Product Name: Recombinant Human NFATC2 Protein
Catalog #: PEH7304



Note

For Research Use Only , Not for Diagnostic Use