


Product Name: Recombinant Human PDE3A Protein
Catalog #: PEH7307



Summary

Name	PDE3A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14432
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSFTSSYAIS AANHVKAKKQ SRPGALAKIS PLSSPCSSPL QGTPASSLVS KISAVQFPES ADTTAKQSLG SHRALTYTQS APDLSPQILT PPVICSSCGR PYSQGNPADE PLERSGVATR TPSRTDDTAQ VTSDYETNNN SDSSDIVQNE DETECLREPL

image



Background

Alternative Names	cGMP-inhibited 3',5'-cyclic phosphodiesterase A (EC 3.1.4.17) (Cyclic GMP-inhibited phosphodiesterase A) (CGI-PDE A)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PDE3A Protein
Catalog #: PEH7307

