


**Product Name: Recombinant Human HSP90AA2P Protein**  
**Catalog #: PEH7309**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | HS902_HUMAN                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                         |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                                                                                            |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Accession #</b>           | Q14568                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Molecular Mass</b>        | 40 kDa                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                          |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                              |
| <b>Sequence</b>              | MPEETQTQDQ PMEEEEVETF AFQAEIAQLM SLIINTFYNS KEIFLRELIS<br>NSSDALDKIW YESLTDPSKL DSGKELHINL IPNKQDQTLT IVDTGIGMTK<br>ADLINNLGTI AKSGTKAFME ALQAGADISM IGQFGVSFYS AYLVAEKVTV<br>ITKHNDDEQY AWESSAGGSF TVRTDTGERM GRGTVKVLHL KEDQTEYLEE<br>QRIKEIVKKH SQLIGYPITL FVEKECDKEV SDDETEEKED KEEEEKEEEK<br>ESKDKPEIED VGSDEEEEEK DGDKKKKKTK EKYIDQEELN KTKPIWTRNP<br>DDITNEEYGE FCKNLTNDWE DHLAVKHFSV EGQLEFRALL FVP |

## image



## Background

|                          |                                                                                                                               |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Heat shock protein HSP 90-alpha A2 (Heat shock 90 kDa protein 1 alpha-like 3) (Heat shock protein HSP 90-alpha A2 pseudogene) |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------|



**Product Name: Recombinant Human HSP90AA2P Protein**  
**Catalog #: PEH7309**

---

**Note**

For Research Use Only , Not for Diagnostic Use