


Product Name: Recombinant Human TTF1 Protein
Catalog #: PEH7314



Summary

Name	TTF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15361
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HTPVSDKKKK KCSIHKERPQ KHSHEIFRDS SLVNEQSQIT RRKKRKKDFQ HLISSPLKKS RICDETANAT STLKKRKKRR YSALEVDEEA GVTVVLVDKE NINNTPKHFR KDVDVVCVDM SIEQKLPRKP KTDKFQVLAK SHAHKSEALH SKVREKKNKK HQRKAASWES QRARDTLPQS ESHQEESWLS VGPGGEITEL PASAHKNKSK KKKKKSSNRE YETLAMPEGS QAGREAGTDM QESQPTVGLD DETPQLLGPT HKKKSKKKKK KSNHQEFEA LAMPEGSQVG SEVGADMQES RPAVGLHGET AGIPAPAYKN KSKKKKKKSN HQEFEAVAMP ESLESAYPEG

image



Background

Alternative Names	Transcription termination factor 1 (TTF-1) (RNA polymerase I termination factor) (Transcription termination factor I) (TTF-I)
--------------------------	---

Product Name: Recombinant Human TTF1 Protein
Catalog #: PEH7314



Note

For Research Use Only , Not for Diagnostic Use