


Product Name: Recombinant Human SYCP1 Protein
Catalog #: PEH7315



Summary

Name	SYCP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15431
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPPRSSSSQV SAVKPQTLGG DSTFFKSFNK CTEDDFEFPF AKTNLSKNGE NIDSDPALQK VNFLPVLEQV GNSDCHYQEG LKDSDLNSE GLSRVYSKLY KEAEKIKKWK VSTEAE LRQK ESKLQENRKI IEAQRKAIQE LQFGNEKVSL KLEEGIQENK DLIKENNATR HLCNLLKETC ARSAEKTCKY EYEREETRQV YMDLNNNIEK MITAFEELRV QAENSRLMH FKLKEDYEKI QHLEQEYKKE INDKEKQVSL LLIQITEKEN KMKDLTFLLE ESRDKVNQLE ETKKLQSEN KQSIEKQHHL TKELEDIKVS LQRSVSTQKA LEEDLQIATK TICQLTEEKE

image



Background

Alternative Names	Synaptonemal complex protein 1 (SCP-1) (Cancer/testis antigen 8) (CT8)
--------------------------	--

Product Name: Recombinant Human SYCP1 Protein
Catalog #: PEH7315



Note

For Research Use Only , Not for Diagnostic Use