


Product Name: Recombinant Human IFI16 Protein
Catalog #: PEH7319



Summary

Name	IF16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16666
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVDATSPAPS TSSTVKTEGA EATPGAQKRK KSTKEKAGPK GSKVSEEQTQ PPSPAGAGMS TAMGRSPSPK TSLAPPNSS STENPKTVAK CQVTPRRNVL QKRPVIVKVL STTKPFYET PEMEKKIMFH ATVATQTQFF HVKVLNTSLK EKFNGKKIII ISDYLEYDSL LEVNEESTVS EAGPNQTFEV PNKIINRAKE TLKIDILHKQ ASGNIVYGVF MLHKKTVNQK TTIYEIQDDR GKMDVVG TGQ CHNIPCEEGL KLQLFCFRLR KKNQMSKLIS EMHSFIQIKK KTNPRNDPK SMKLPQEQRQ LPYPSEASTT FPESHLRTPQ MPPTTPSSSF

image



Background

Alternative Names Gamma-interferon-inducible protein 16 (I γ I-16) (Interferon-inducible myeloid differentiation transcriptional activator)

Product Name: Recombinant Human IFI16 Protein
Catalog #: PEH7319



Note

For Research Use Only , Not for Diagnostic Use