


Product Name: Recombinant Human CYP1B1 Protein
Catalog #: PEH7320



Summary

Name	CP1B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16678
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRELLSHNEE FGRTVGAGSL VDVMPWLQYF PNPVRTVFRE FEQLNRNFSN FILDKFLRHC ESLRPGAAPR DMMDAFILSA EKKAAGDSHG GGARLDLENV PATITDIFGA SQDTLSTALQ WLLLLFTRYP DVQTRVQAEI DQVVGRDRLP CMGDQPNLPY VLAFLYEAMR FSSFVPVTIP HATTANTSVL GYHIPKDTVV FVNQWSVNHD PLKWPNPENF DPARFLDKDG LINKDLTSRV MIFSVGKRRC

image



Background

Alternative Names Cytochrome P450 1B1 (EC 1.14.14.1) (CYP1B1) (Hydroperoxy
icosatetraenoate dehydratase) (EC 4.2.1.152)

Note

Product Name: Recombinant Human CYP1B1 Protein
Catalog #: PEH7320



For Research Use Only , Not for Diagnostic Use