


Product Name: Recombinant Human SV2A Protein
Catalog #: PEH7321



Summary

Name	SV2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L0J3
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FIRGAKDIAK EVKKHAAKKV VKGLDRVQDE YSRYSRFE EEDDDDFPA PSDGYRGEQ TQDEEEGGAS SDATEGHDED DEIYEGEYQG IPRAESGGKG ERMADGAPLA GVRGGLSDGE GPPGGRGEAQ RRKEREELAQ QYEAILRECG

image



Background

Alternative Names Synaptic vesicle glycoprotein 2A

Note

For Research Use Only , Not for Diagnostic Use