

Product Name: Recombinant Human ASRGL1 Protein
Catalog #: PEH7322



Summary

Name	ASGL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L266
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNPIVVHGG GAGPISKDRK ERVHQGMVRA ATVGYGILRE GGSVDAVEG AVVALEDDPE FNAGCGSVLN TNGEVEMDAS IMDGKDLAS AVSAVQCIAN PIKLARLVME KTPHCFLTDQ GAAQFAAAMG VPEIPGEKLV TERNKKRLEK EKHEKGAQKT DCQKNLGTVG AVALDCKGNV AYATSTGGIV NKMVGRVGD PCLGAGGYAD NDIGAVSTTG HGESILKVNL ARLTLFHIEQ GKTVEEAADL SLGYMKSRVK GLGLIVVSK TGDWVAKWTS TSMPWAAAKD GKLHFGIDPD DTTITDLP

image



SDS



WB

Background

Alternative Names Isoaspartyl peptidase/L-asparaginase (EC 3.4.19.5) (EC 3.5.1.1)
(Asparaginase-like protein 1) (Beta-aspartyl-peptidase) (Isoaspartyl

Product Name: Recombinant Human ASRGL1 Protein
Catalog #: PEH7322



dipeptidase) (L-asparagine amidohydrolase) [Cleaved into: Isoaspartyl
peptidase/L-asparaginase alpha chain; Isoaspartyl peptidase/L-asparaginase
beta chain]

Note

For Research Use Only , Not for Diagnostic Use