


Product Name: Recombinant Human PIAS4 Protein
Catalog #: PEH7325



Summary

Name	PIAS4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N2W9
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSELLQRLKT IGVKHPSELCK ALVKEKLRLD PDSEIATTGV RVSLICPLVK MRLSVPCRAE TCAHLQCFDA VFYLMNEKK PTWMCPVCDK PAPYDQLIID GLLSKILSEC EDADEIEYLV DGSWCPIRAE KERSCSPQGA ILVLGPSDAN GLLPAPSVNG SGALGSTGGG GPVGSMENGK PGADVVDLTL DSSSSSEDEE EEEEEEEEDED EEGPRPKRRC

image



Background

Alternative Names	E3 SUMO-protein ligase PIAS4 (EC 2.3.2.27) (PIASy) (Protein inhibitor of activated STAT protein 4) (Protein inhibitor of activated STAT protein gamma) (PIAS-gamma) (RING-type E3 ubiquitin transferase PIAS4)
--------------------------	--

Product Name: Recombinant Human PIAS4 Protein
Catalog #: PEH7325



Note

For Research Use Only , Not for Diagnostic Use