


Product Name: Recombinant Human SYT2 Protein
Catalog #: PEH7326



Summary

Name	SYT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N9I0
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETGLTEGEGE GEEKEPENL GKLFSLDYD FQANQLTVGV LQAAELPALD MGGTSDPYVK VFLLPDKKKK YETKVHRKTL NPAFNETFTF KVPYQELGGK TLVMAIYDFD RFSKHDIIGE VKVPMNTVDL GQPIEEWRDL QGGEKEEPEK LGDICTSLRY VPTAGKLTVC ILEAKNLKKM DVGGLSDPYV KIHLMQNGKR LKKKKTTVKK KTLNPYFNES FSFEIPFEQI QKVQVVVTVL DYDKLGKNEA IGKIFVGSNA

image



Background

Alternative Names Synaptotagmin-2 (Synaptotagmin II) (SytII)

Product Name: Recombinant Human SYT2 Protein
Catalog #: PEH7326



Note

For Research Use Only , Not for Diagnostic Use