


Product Name: Recombinant Human IPO4 Protein
Catalog #: PEH7329



Summary

Name	IPO4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEX9
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVEEEDDSEI SGYSVENAFF DEKEDTCAAV GEISVNTSVA FLPYMESVFE EVFKLLECPH LNVKAAHEA LGQFCCALHK ACQSCPSEPN TAALQAALAR VVPSYMQAVN RERERQVMA VLEALTGVLR SCGTLTLKPP GRLAELCGVL KAVLQRKTAC QDTDEEEEE

image



Background

Alternative Names Importin-4 (Imp4) (Importin-4b) (Imp4b) (Ran-binding protein 4) (RanBP4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IPO4 Protein
Catalog #: PEH7329

