


Product Name: Recombinant Human ABCC2 Protein
Catalog #: PEH7331



Summary

Name	MRP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q92887
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EATVHDGSEE EDDDYGLISS VEEIPEDAAS ITMRRENSFR RTLSRSSRSN GRHLKSLRNS LKTRNVNSLK EDEELVKGQK LIKKEFIETG

image



Background

Alternative Names	Canalicular multispecific organic anion transporter 1 (ATP-binding cassette sub-family C member 2) (Canalicular multidrug resistance protein) (Multidrug resistance-associated protein 2) (EC 7.6.2.2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use