


Product Name: Recombinant Human IPO9 Protein
Catalog #: PEH7337



Summary

Name	IPO9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96P70
Host	E.coli
Species	Human
Molecular Mass	72 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISKAGRELGE NLDQILRAIL SKMQQAETLS VMQSLIMVFA HLVHTQLEPL LEFLCSLPGP TGKPALEFVM AEWTSRQHLF YGQYEGKVSS VALCKLLQHG INADDKRLQD IRVKGEEIYS MDEGIRTRSK SAKNPERWTN IPLLVKILKL IINELSNVME ANAARQATPA EWSQDDSNM WEDQEEEEEE EEDGLAGQLL SDILATSKYE EDYYEDDEED DPDALKDPLY QIDLQAYLTD FLCQFAQQPC YIMFSGHLND

image



Background

Alternative Names Importin-9 (Imp9) (Ran-binding protein 9) (RanBP9)

Product Name: Recombinant Human IPO9 Protein
Catalog #: PEH7337



Note

For Research Use Only , Not for Diagnostic Use