


Product Name: Recombinant Human RNF5 Protein
Catalog #: PEH7339



Summary

Name	RNF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q99942
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEGPNRERGG AGATFECNIC LETAREAVVS VCGHLYCWPC LHWLETRPE RQECPVCKAG ISREKVVPLY GRGSQKPQDP RLKTPRPQGG QRPAPESRGG

image



Background

Alternative Names E3 ubiquitin-protein ligase RNF5 (EC 2.3.2.27) (Protein G16) (RING finger protein 5) (RING-type E3 ubiquitin transferase RNF5) (Ram1 homolog) (HsRma1)

Note

For Research Use Only , Not for Diagnostic Use