


Product Name: Recombinant Human TEP1 Protein
Catalog #: PEH7340



Summary

Name	TEP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99973
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LNFDINLENP SRTLISITQA KPESESSFLC ASSDGILWNL AKCSPEGEWT TGNMWQKKAN TPETQTPGTD PSTCRESNAS MDSNASMDSE PTPHLKTRQR RKIHSGSVTA LHVLPPELLVT

image



Background

Alternative Names	Telomerase protein component 1 (Telomerase-associated protein 1) (Telomerase protein 1) (p240) (p80 telomerase homolog)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use