


Product Name: Recombinant Human MAPKAP1 Protein
Catalog #: PEH7341



Summary

Name	SIN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BPZ7
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LIDHDVDLEK IHPSPMPGDS GSEIQGSNGE TQGYVYAQSV DITSSWDFGI RRRSNTAQRL ERLRKERQNN IKCKNIQWKE RNSKQSAQEL KSLFEKSLK EKPPISGKQS ILSVRLEQCP LQLNPFNEY SKFDGKGHVG TTATKKIDVY LPLHSSQDRL LPMTVVVTMAS ARVQDLIGLI CWQYTSEGRE PKLNDNVSAY CLHIAEDDGE VDTDFPPLDS NEPIHKFGFS

image



Background

Alternative Names	Target of rapamycin complex 2 subunit MAPKAP1 (TORC2 subunit MAPKAP1) (Mitogen-activated protein kinase 2-associated protein 1) (Stress-activated map kinase-interacting protein 1) (SAPK-interacting protein 1) (mSIN1)
--------------------------	--

Product Name: Recombinant Human MAPKAP1 Protein
Catalog #: PEH7341



Note

For Research Use Only , Not for Diagnostic Use