


Product Name: Recombinant Human HDAC6 Protein
Catalog #: PEH7349



Summary

Name	HDAC6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBN7
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAACTRSLLG DPPLLTLPR PPLSGALASI TETIQVHRRY WRSLRVMKVE DREGPSSSKL VTKKAPQPAK PRLAERMTR EKKVLEAGMG KVTSASFGEE STPGQTNSSET AVVALTQDQP SEAATGGATL AQTISEAAIG GAMLGQTTSE EAVGGATPDQ TTSEETVGGGA ILDQTTSEDA VGGATLGQTT SEEAVGGATL AQTTEAAME GATLDQTTSE EAPGGTELIQ TPLASSTDHQ TPPTSPVQGT TPQISPSTLI GSLRTLELGS ESQGASESQA PGEENLLGEA AGGQDMADSM LMQGSRGLTD QAIFYAVTPL

image



Background

Alternative Names Histone deacetylase 6 (HD6) (EC 3.5.1.98) (Tubulin-lysine deacetylase HDAC6) (EC 3.5.1.-)

Product Name: Recombinant Human HDAC6 Protein
Catalog #: PEH7349



Note

For Research Use Only , Not for Diagnostic Use