


Product Name: Recombinant Human ASIC3 Protein
Catalog #: PEH7351



Summary

Name	ASIC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHC3
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLDCFRGQP CGPENFTTIF TRMGKCYTFN SGADGAELLT TTRGGMGNGL DIMLDVQEE YLPVWRDNEE TPFEVGIRVQ IHSQEEPPII DQLGLGVSPG

image



Background

Alternative Names	Acid-sensing ion channel 3 (ASIC3) (hASIC3) (Amiloride-sensitive cation channel 3) (Neuronal amiloride-sensitive cation channel 3) (Testis sodium channel 1) (hTNaC1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use