


Product Name: Recombinant Human CROT Protein
Catalog #: PEH7353



Summary

Name	OCTC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9UKG9
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERTFQYQDSL PSLPVPSLEE SLKKYLESVK PFANQEEYKK TEEIVQKFQS GIGEKLHQKL LERAKGKRNW LEEWWLNVAY LDVRIPSQLN VNFAGPAAHF EHYWPPKEGT QLERSITLW HNLNYWQLLR KEKVPVHKVG NTPLDMNQFR MLFSTCKVPG ITRDSIMNYF RTESEGRSPN HIVVLCRGRA FVFDVIHEGC LVTPPELLRQ LTYIHKKCHS EPDGPZIAAL TSEERTRWAK AREYLIGLDP

image



Background

Alternative Names Peroxisomal carnitine O-octanoyltransferase (COT) (EC 2.3.1.137)

Note

Product Name: Recombinant Human CROT Protein
Catalog #: PEH7353



For Research Use Only , Not for Diagnostic Use