


Summary

Name	PDE10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y233
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEVSRYQDTN MQGVVYELNS YIEQRLDTGG DNQLLYELS SIIKIATKAD GFALYFLGEC NNSLCIFTPP GIKEGKPRLI PAGPITQGTT VSAYVAKSRK TLLVEDILGD ERFPRGTGLE SGTRIQSVLC LPIVTAIGDL IGILELYRHW GKEAFCLSHQ EVATANLAWA SVAIHQVQVC RGLAKQTELN DFLLDVSKTY

image



Background

Alternative Names cAMP and cAMP-inhibited cGMP 3',5'-cyclic phosphodiesterase 10A (EC 3.1.4.17) (EC 3.1.4.35)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PDE10A Protein
Catalog #: PEH7357

