


Product Name: Recombinant Human SCN10A Protein
Catalog #: PEH7358



Summary

Name	SCNAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5Y9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEQNQATTDE IEAKEKKFQE ALEMLRKEQE VLAALGIDTT SLHSHNGSPL TSKNASERRH RIKPRVSEGS TEDNKSPRS PYNQRRMSFL GLASGKRRAS HGVSFHFRRSP GRDISLPEGV TDDGVFPGDH ESHRGSLLLG

image



Background

Alternative Names	Sodium channel protein type 10 subunit alpha (Peripheral nerve sodium channel 3) (PN3) (hPN3) (Sodium channel protein type X subunit alpha) (Voltage-gated sodium channel subunit alpha Nav1.8)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SCN10A Protein
Catalog #: PEH7358

