


Product Name: Recombinant Human MYH4 Protein
Catalog #: PEH7359



Summary

Name	MYH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y623
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KMEGDLNEME IQLNHANRQA AEALRNLRNT QGILKDTQLH LDDAIRGQDD LKEQLAMVER RANLMQAEVE ELRASLERTE RGRKMAEQEL LDASERVQLL HTQNTSLINT KKKLETDISQ IQGEMEDIVQ EARNAEKAK KAITDAAMMA EELKKEQDTS AHLERMKKNM EQTVKDLQLR LDEAEQLALK GGKKQIQKLE ARVRELESEV ESEQKHNEVA VKGLRKHERR VKELTYQTEE DRKNILRLQD LVDKLQTKVK AYKRQAEAE EQSNVNLAKF RKLQHELEEA KERADIAESQ VNKL RVKSRE

image



Background

Alternative Names	Myosin-4 (Myosin heavy chain 2b) (MyHC-2b) (Myosin heavy chain 4) (Myosin heavy chain IIb) (MyHC-IIb) (Myosin heavy chain, skeletal muscle,
--------------------------	---

Product Name: Recombinant Human MYH4 Protein
Catalog #: PEH7359



fetal)

Note

For Research Use Only , Not for Diagnostic Use