


Product Name: Recombinant Human PGP Protein
Catalog #: PEH7362



Summary

Name	PGP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NDG6
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAAAEAGGDD ARCVRLSAER AQALLADVDT LLFDCDGV LW RGETAVPGAP EALRALRARG KRLGFITNNS SKTRAAYA EK LRRLGFGGPA GPGASLEVFG TAYCTALYLR QRLAGAPAPK AYVLGSPALA AELEAVGVAS VGVGPEPLQG EGPGDWLHAP LEPDVRAVVV GFDPHFSYMK LTKALRYLQQ PGCLLVGTNM DNRLPLENGR FIAGTGCLVR AVEMAAQRQA DIIGKPSRFI FDCVSQEYGI NPERTVMVGD RLDTDILLGA TCGLKTILTL TGVSTLGDVK NNQESDCVSK KKMVPDFYVD SIADLLPALQ G

image



Background

Alternative Names	Glycerol-3-phosphate phosphatase (G3PP) (EC 3.1.3.21) (Aspartate-based ubiquitous Mg(2+)-dependent phosphatase) (AUM) (EC 3.1.3.48)
--------------------------	---

Product Name: Recombinant Human PGP Protein
Catalog #: PEH7362



(Phosphoglycolate phosphatase) (PGP)

Note

For Research Use Only , Not for Diagnostic Use