


Product Name: Recombinant Human NCF1B Protein
Catalog #: PEH7363



Summary

Name	NCF1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NI72
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGDTFIRHIA LLGFEKRFVP SQHYVRYMFL VKWQDLSEKV VYRRFTEIYE FHKTLKEMFP IEAGAINPEN RIIPHLPAK WFDGQRAAEN HQGTLTEYCG TLMSLPTKIS RCPHLLDFFK VRPDDLKLPT DNQTKKPETY LMPKDGKSTA TDITGPILQ TYRAIANYEK TSGSEMALST GDVVEVVEKS ESGWWFCQMK AKRGWIPASF LEPLDSPDET EDPEPNYAGE PYVAIKAYTA VEGDEVSLLE GEAVEVIHKL LDGWVWIRKD DVTGYFPSMY LQKSGQDVSQ AQRQIKRGAP PRRSSIRNVH SIHQRSRKRL SQDAYRRNSV RFLQQRRRQA RPGPQSPGSP LEERQTQRS KPQPAVPPRP SADLILNRCS ESTKRKLASA V

image



Background

Alternative Names	Putative neutrophil cytosol factor 1B (NCF-1B) (Putative SH3 and PX domain-
--------------------------	---

Product Name: Recombinant Human NCF1B Protein
Catalog #: PEH7363



containing protein 1B)

Note

For Research Use Only , Not for Diagnostic Use