


Product Name: Recombinant Human SLC28A1 Protein
Catalog #: PEH7365



Summary

Name	S28A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00337
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SISLTPVAKG LENMGADFLE SLEEGQLPRS DLSPAEIRSS WSEAAPKPFS RWRNLQPALR ARSFCREHMQ

image



Background

Alternative Names	Sodium/nucleoside cotransporter 1 (Concentrative nucleoside transporter 1) (CNT 1) (hCNT1) (Na(+))/nucleoside cotransporter 1) (Sodium-coupled nucleoside transporter 1) (Solute carrier family 28 member 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use