


**Product Name: Recombinant Human IPO5 Protein**  
**Catalog #: PEH7366**



## Summary

<b>Name</b>	IPO5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O00410
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	35 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RLLRGKTIEC ISLIGLAVGK EKFMQDASDV MQLLLKTQTD FNDMEDDDPQ ISYMISAWAR MCKILGKEFQ QYLPVVMGPL MKTASIKPEV ALLDTQDMEN MSDDDGWEFV NLGDQQSFGI KTAGLEEKST ACQMLVCYAK ELKEGFVEYT EQVVKLMVPL LKFYFHDGVR VAAAESMPLL LECARVRGPE YLTQMWHFMC DALIKAIGTE PDSVLSEIM HSFACIEVM GDGCLNNEHF EELGGILKAK LEEHFKNQEL RQVKRQDEDY DEQVEESLQD

## image



## Background

<b>Alternative Names</b>	Importin-5 (Imp5) (Importin subunit beta-3) (Karyopherin beta-3) (Ran-binding protein 5) (RanBP5)
--------------------------	---

**Product Name: Recombinant Human IPO5 Protein**  
**Catalog #: PEH7366**



---

**Note**

For Research Use Only , Not for Diagnostic Use