


Product Name: Recombinant Human VWA5A Protein
Catalog #: PEH7369



Summary

Name	VMA5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00534
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HREPVLKSI SVSVNIYEFV AGVSATLNYE NEEKVPLEAF VFVPMDEDSA VYSFEALVDG KKIVAELQDK MKARTNYEKA ISQGHQAFLL EGDSSSRDVF SCNVGNLQPG SKAAVTLKYV QELPLEADGA LRFVLPVAVLN PRYQFSGSSK DSCLVNKTPI VPVEDLPYTL SMVATIDSQH GIEKVQSNCP LSPTEYLGED K TSAQVSLAA GHKFD RDVEL LIYYNEVHTP SVVLEMGMMPN MKPGHLMGDP

image



Background

Alternative Names	von Willebrand factor A domain-containing protein 5A (Breast cancer suppressor candidate 1) (BCSC-1) (Loss of heterozygosity 11 chromosomal region 2 gene A protein)
--------------------------	--

Product Name: Recombinant Human VWA5A Protein
Catalog #: PEH7369



Note

For Research Use Only , Not for Diagnostic Use