

Product Name: Recombinant Human TRAFD1 Protein
Catalog #: PEH7372



Summary

Name	TRAD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14545
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLCDNCKKEI PVFNFTIHEI HCQRNIGMCP TCKEFPKSD METHMAAEHC QVTCKCNKKL EKRLKKHEE TECPLRLAVC QHCDLELSIL KLKEHEDYCG ARTELCGNCG RNVLVKDLKT HPEVCGREGE EKRNEVAIPP NAYDESWGQD GIWIASQLLR QIEALDPPMR LPRRPLRAFE SDVFHNRTTN QRNITAQVSI QNNLFEEQER QERNRGQQPP KEGGEESANL DFMLALSLQN EGQASSVAEQ DFWRAVCEAD QSHGGPRSLS DIKGADEIM LPCEFCEELY PEELLIDHQT SCNPSRALPS LNTGSSSPRG VEEDVIFQN

image



Background

Alternative Names TRAF-type zinc finger domain-containing protein 1 (Protein FLN29)

Product Name: Recombinant Human TRAFD1 Protein
Catalog #: PEH7372



Note

For Research Use Only , Not for Diagnostic Use