


Product Name: Recombinant Human AP3D1 Protein
Catalog #: PEH7373



Summary

Name	AP3D1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14617
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPLKVPGLPM SDQYVKLEEE RRHRQKLEKD KRRKKRKEKE KKGKRRHSSL PTESDEDIAP AQQVDIVTEE MPENALPSDE DDKDPNDPYR ALDIDLKPL ADSEKLPIQK HRNTETSKSP

image



Background

Alternative Names AP-3 complex subunit delta-1 (AP-3 complex subunit delta) (Adaptor-related protein complex 3 subunit delta-1) (Delta-adaptin)

Note

For Research Use Only , Not for Diagnostic Use