


Product Name: Recombinant Human ERC2 Protein
Catalog #: PEH7378



Summary

Name	ERC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15083
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEKERIIERL KEQRERDDRE RLEEIESFRK ENKDLKEKVN ALQAELTEKE SSLIDLKEHA SSLASAGLKR DSKLKSLEIA IEQKKEECKS LEAQLKKAHN IEDDSRMNPE FADQIKQLDK EASYRDECG KAQAEVDRLL EILKEVENEK NDKDKKIAEL ESLTLRHMKD QNKKVANLKH NQQLKKKNA QLLEEVRRRE DSMADNSQHL QIEELMNALE KTRQELDATK ARLASTQQSL AEKEAHLANL RIERRKQLEE ILEMKQEALL AAISEKDANI ALLELSASKK KKTQEVMAL KREKDRLVHQ LKQQTQNRMK LMADNYDDDH HHYHHHHHHH HHRSPGRSQH SNHRPSPDQD

image



Background

Alternative Names ERC protein 2

Product Name: Recombinant Human ERC2 Protein
Catalog #: PEH7378



Note

For Research Use Only , Not for Diagnostic Use