


Summary

Name	SMAD7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O15105
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WCVVAYWEEK TRVGRLYCVQ EPSLDIFYDL PQGNGFCLGQ LNSDNKSQLV QKVRISKIGCG IQLTREVDGV WVYNRSSYPI FIKSATLDNP DSRTLLVHKV FPGFSIKAFD YEKAYSLQRP NDHEFMQQPW TGFTVQISFV KGWGQCYTRQ FISSPCWLE

image



Background

Alternative Names	Mothers against decapentaplegic homolog 7 (MAD homolog 7) (Mothers against DPP homolog 7) (Mothers against decapentaplegic homolog 8) (MAD homolog 8) (Mothers against DPP homolog 8) (SMAD family member 7) (SMAD 7) (Smad7) (hSMAD7)
--------------------------	--

Product Name: Recombinant Human SMAD7 Protein
Catalog #: PEH7380



Note

For Research Use Only , Not for Diagnostic Use