


Summary

Name	P4HA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15460
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAAKALMRLQ DTYRLDPGTI SRGELPGTKY QAMLSVDDCF GMGRSAYNEG DYYHTVLWME QVLKQLDAGE EATTTKSQVL DYLSYAVFQL GDLHRALELT RRLSLDPSH ERAGGNLRYF EQLLEEEREK TLTNQTEAEL ATPEGIYERP VDYLPERDVY

image



Background

Alternative Names	Prolyl 4-hydroxylase subunit alpha-2 (4-PH alpha-2) (EC 1.14.11.2) (Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human P4HA2 Protein
Catalog #: PEH7383

