

Product Name: Recombinant Human SART1 Protein
Catalog #: PEH7387



Summary

Name	SNUT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43290
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEAAGTTAAA GTGGATEQPP RHREHKHKH RSGSGSGSGG ERRKRSRERG GERGSGRRGA EAEARSSTHG RERSQAEPSE RRVKREKRDD GYEAAASSKT SSGDASSLSI EETNKLRAKL GLKPLEVNAI KKEAGTKEEP VTADVNPMA LRQREELREK LAAAKEKRLL NQKLGIKTL GEDDPWLDDT AAWIERSRQL QKEKDLAEKR AKLLEEMDQE FGVSTLVEEE FGQRRQDLYS ARDLQGLTVE HAIDSFREGE TMILTLKDKG VLQEEEDVLV NVNLVDKERA EKNVELRKKK PDYLPYAEDE SVDDLAQQKP RSILSKYDEE LEGERPHSFR LEQGGTADGL RERELEEIRA KLRLQAQSL

image



Background

Alternative Names	U4/U6.U5 tri-snRNP-associated protein 1 (SNU66 homolog) (hSnu66)
--------------------------	--

Product Name: Recombinant Human SART1 Protein
Catalog #: PEH7387



(Squamous cell carcinoma antigen recognized by T-cells 1) (SART-1)
(hSART-1) (U4/U6.U5 tri-snRNP-associated 110 kDa protein) (allergen Hom s
1)

Note

For Research Use Only , Not for Diagnostic Use