


Product Name: Recombinant Human MTSS1 Protein
Catalog #: PEH7388



Summary

Name	MTSS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43312
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALGGLFQTII SDMKGSYPVW EDFINKAGKL QSQLRTTVVA AAFLDAFQK VADMATNTRG GTREIGSALT RMCMRHRSIE AKLRQFSSAL IDCLINPLQE QMEEWKKVAN QLDKDHAKY KKRQEIKKK SSDTLKLQKK AKKGRGDIQP QLDSALQDVN DKYLLLEETE KQAVRKALIE ERGRFCTFIS MLRPVIEEEI SMLGEITHLQ TISEDLKS LT MDPHKLPS SS EQVILD LKGS

image



Background

Alternative Names Protein MTSS 1 (Metastasis suppressor YGL-1) (Metastasis suppressor protein 1) (Missing in metastasis protein)

Note

Product Name: Recombinant Human MTSS1 Protein
Catalog #: PEH7388



For Research Use Only , Not for Diagnostic Use