


Product Name: Recombinant Human EIF4G3 Protein
Catalog #: PEH7392



Summary

Name	IF4G3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43432
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQELFRKQRS ILNKLTPQMF NQLMKQVSGL TVDTEERLKG VIDLVFEKAI DEPSFSVAYA NMCRCVLTK VPMADKPGNT VNFRKLLNR CQKEFEKDKA DDDVFEKKQK ELEAASAPEE RTRLHDELEE AKDKARRRSI GNIKFIGELF KLKMLTEAIM HDCVVKLLKN HDEESLECLC RLLTTIGKDL DFEKAKPRMD QYFNQMEKIV KERKTSSRIR FMLQDVIDLR LCNWVSRRAD QGPKTIEQIH

image



Background

Alternative Names Eukaryotic translation initiation factor 4 gamma 3 (eIF-4-gamma 3) (eIF-4G 3) (eIF4G 3) (eIF-4-gamma II) (eIF4GII)

Note

Product Name: Recombinant Human EIF4G3 Protein
Catalog #: PEH7392



For Research Use Only , Not for Diagnostic Use