


Product Name: Recombinant Human LAMTOR5 Protein
Catalog #: PEH7393



Summary

Name	LTOR5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43504
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEATLEQHLE DTMKNPSIVG VLCTDSQGLN LGCRGTLSD E HAGVISVLAQ QAAKLTSDPT DIPVVCLESD NGNIMI QKHD GITVAVHKMA S

image



Background

Alternative Names	Ragulator complex protein LAMTOR5 (Hepatitis B virus X-interacting protein) (HBV X-interacting protein) (HBX-interacting protein) (Late endosomal/lysosomal adaptor and MAPK and MTOR activator 5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use